










Accu-Chek® Mobile – Proper Testing Instructions

Incorrect handling may lead to incorrect test results which can cause the wrong therapy recommendation to be made and so result in serious adverse health effects. Therefore, follow these instructions carefully to perform a blood glucose test correctly.

Do	Do not	Why
<p>1</p>  <p>Wash your hands with warm water and soap and rinse well. Dry your hands thoroughly with a clean towel before obtaining blood.</p>	 <p>Do not test when your fingers are unwashed, soiled or sticky (e. g. due to food or drink residue).</p>	<p>1</p> <p>Unwashed, soiled or sticky fingers can lead to falsely elevated blood glucose results and delivery of an inappropriately high insulin amount.</p>
<p>2</p>  <p>Always use a small, fresh and well-formed blood drop.</p>	 <p>Do not test with a smeared or spread blood drop. Do not smear the blood drop onto or across the test area.</p>	<p>2</p> <p>Testing with a spread blood drop or wiping the blood drop on the test area increases the risk of contamination of the testing site. This can lead to falsely elevated blood glucose results and delivery of an inappropriately high insulin amount.</p>
<p>3</p>  <p>Touch the test area gently and steadily when applying blood. Apply the blood drop to the centre of the test area.</p>	 <p>Do not press the test area onto the finger (especially if shaking or trembling).</p>	<p>3</p> <p>Pressing the test area onto the finger can lead to falsely elevated blood glucose results and delivery of an inappropriately high insulin amount.</p>
<p>4</p> <p>Apply the blood drop to the test area immediately after you have obtained the blood drop. Remove your finger from the test area when the beep tone sounds and <i>Test in progress</i> is displayed.</p> <p>Note: You can only hear beep tones when they are turned on. See User's Manual for tone settings.</p> 	<p>Do not wait before applying the blood drop to the test area. Do not hold your finger on the test area after the beep tone sounds and <i>Test in progress</i> is displayed.</p>	<p>4</p> <p>Not testing promptly can lead to evaporation of the blood drop. This can lead to falsely elevated blood glucose results and delivery of an inappropriately high insulin amount. The beep tone helps you to obtain a reliable test result.</p>

www.accu-chek.com.au | Accu-Chek Enquiry Line: 1800 251 816

ACCUCHEK®

For people with diabetes. Use only as directed. Consult your healthcare professional for advice. The Accu-Chek Mobile test cassettes are PBS-listed for insulin-requiring people with diabetes. ACCU-CHEK and ACCU-CHEK MOBILE are trademarks of Roche. © 2014 Roche Diagnostics. All other trademarks are the property of their respective owners.

Roche Diagnostics Australia Pty Ltd
31 Victoria Ave, Castle Hill, NSW, 2154 ABN 29 003 001 205